

## Gypsum Analysis

### 1. **Scope:**

This document provides a standardized procedure for the analysis of natural gypsum or gypsite. This method is not suitable for the analysis of by-product gypsums that contain free sulfuric acid and/or other sulfates. This method may be used to determine sulfate in other products such as potassium sulfate or magnesium sulfate.

### 2. **Principle:**

The ground sample is dissolved in dilute HCl and the sulfate is precipitated as BaSO<sub>4</sub>. The percent gypsum is calculated from the mass of BaSO<sub>4</sub>. All analyses are reported based on the material as received and not on a dry weight basis.

### 3. **Safety:**

Read the SDS for all materials before use.

### 4. **Apparatus and Equipment:**

Equivalents are acceptable

- 4.1. Aluminum moisture dishes with covers (VWR cat. # 25420-072)
- 4.2. Drying Ovens – 45°C and 105°C
- 4.3. 400mL beakers
- 4.4. Hot plate or steam bath
- 4.5. Whatman binder free glass microfiber filter, Type (GF/F), Particle retention 0.7 µm (Whatman Cat. #1825-024)
- 4.6. Porcelain Gooch Crucibles (Coors Cat. # 60148 or 60151)
- 4.7. Vacuum set up
- 4.8. Analytical balance
- 4.9. Desiccator

### 5. **Reagents and Supplies:**

- 5.1. Hydrochloric acid, conc.
- 5.2. Barium chloride

### 6. **Free Water Determination:**

- 6.1. Prepare samples as described in RA-SP-SMPL-PREP.
- 6.2. Mix sample thoroughly and weigh 20.00g into an aluminum dish. Dry at 45°C for a minimum of 8 hours or overnight. Keep covered and store in a desiccator to cool.
- 6.3. Re-weigh the cooled sample and calculate the weight loss as the moisture (or “free water”) in the material.

7. **BaSO<sub>4</sub> Determination:**

- 7.1. Prepare 10% barium chloride solution by adding 100g BaCl<sub>2</sub> to a 1000mL volumetric flask. Add water to dissolve then fill to the mark. Filter through Whatman No. 42 paper.
- 7.2. Weigh 0.5 ±0.001g of the ground sample into a 400mL beaker.
- 7.3. Add 10mL HCl and then 50mL H<sub>2</sub>O. Cover with a watch glass and gently boil for ~10 minutes.
- 7.4. Filter boiling solution with vacuum through a glass microfiber filter in a Gooch crucible. Rinse beaker with hot water and pour the rinse onto the filter. Repeat as needed to be certain that all material is washed onto the filter.
- 7.5. Thoroughly clean beaker with a brush and rinse with DI water. Transfer filtrate back to cleaned beaker (volume should be approximately 250mL).
- 7.6. Cover with a watch glass and bring to a boil. Add 15mL of 10% BaCl<sub>2</sub> dropwise until solution is well-clouded (BaSO<sub>4</sub> forms as the precipitate). The remainder may be added more rapidly. Keep on a hot steam bath or keep warm at low temperature for at least one hour. Allow to cool for several hours.
- 7.7. Dry an additional clean Gooch crucible with a glass microfiber filter at ~105°C for at least one hour. Cool in a desiccator and record the weight. Filter cooled solution with vacuum through the crucible. Wash filter at least 5 times with room temperature DI water.
- 7.8. Dry crucible and contents for 2 hours at ~105°C. Cool in a desiccator and weigh to obtain the mass of BaSO<sub>4</sub>.

8. **Equations:**

$$\% \text{ Gypsum} = \frac{\text{mass of BaSO}_4 \times 0.7377}{\text{mass of sample}} \times 100$$

$$\% \text{ Sulfur} = \frac{\text{mass of BaSO}_4 \times 0.1374}{\text{mass of sample}} \times 100$$

$$\% \text{ Gypsum} = \frac{\% \text{ Sulfur}}{0.1863}$$

$$\% \text{ Gypsum} = \frac{\% \text{ Calcium}}{0.2328}$$

Correction for loss of moisture during grinding:

$$\% \text{ Gypsum, as received} = \frac{100 - \% \text{ moisture (unground)}}{100 - \% \text{ moisture (ground)}} \times \% \text{ Gypsum (ground)}$$

Note:

$$0.7377 = \frac{\text{Molecular weight Gypsum}}{\text{Molecular weight BaSO}_4} = \frac{172.17}{233.38}$$

$$0.1374 = \frac{\text{Molecular weight Sulfur}}{\text{Molecular weight BaSO}_4} = \frac{32.07}{233.38}$$

$$0.1863 = \frac{\text{Molecular weight Sulfur}}{\text{Molecular weight Gypsum}} = \frac{32.07}{172.17}$$

$$0.2328 = \frac{\text{Molecular weight Calcium}}{\text{Molecular weight Gypsum}} = \frac{40.08}{172.17}$$

Gypsum equivalent is calculated both using the mass of BaSO<sub>4</sub> and the percent calcium from mineral analysis and both values are reported.

Results and reruns are reported per client's request.

**Approvals:**

**Revised By:**

*Stacy Aylesworth*  
Stacy Aylesworth  
Senior Environmental Scientist

6/18/20  
Date

**Approved By:**

*Maryam Khosravifard*  
Maryam Khosravifard  
Environmental Program Manager I

6/18/20  
Date

**Approved By:**

*Sarva Balachandra*  
Sarva Balachandra  
Quality Assurance Officer

6/28/20  
Date

