

## **Moisture / Free Water Analysis in Fertilizer**

### **1. Scope:**

This document provides a standardized procedure for determining the moisture or free water in fertilizer samples.

### **2. Principle:**

The ground sample is dried in an oven and the difference between the initial weight and the dried weight is used to calculate % moisture or free water. Samples are prepared as described in RA-SP-SMPL-PREP. Gypsum samples should not be pre-dried in an oven prior to analysis.

### **3. Safety:**

Care should be taken when removing hot items from ovens.

### **4. Equipment:**

- 4.1. Analytical balance
- 4.2. Drying Ovens – 45°C and 105°C
- 4.3. Vacuum oven – 65°C
- 4.4. Desiccator
- 4.5. Aluminum moisture dishes with covers (VWR cat. # 25420-072)

### **5. Moisture Determination in Solid Samples:**

- 5.1. Weigh the empty aluminum dish and record the weight.
- 5.2. Mix sample thoroughly and weigh  $2.0\text{g} \pm 0.1\text{g}$  into the aluminum dish (record the weight).
- 5.3. Dry at 105°C for ~3 hours. Keep covered and store in a desiccator to cool.
- 5.4. Re-weigh the cooled sample and calculate the weight loss as the moisture.

## 6. Moisture Determination in Liquid Samples:

- 6.1. Add ~10.0g sand to an aluminum dish and record the weight.
- 6.2. Mix sample thoroughly and weigh  $2.0\text{g} \pm 0.1\text{g}$  into the aluminum dish (record the weight).
- 6.3. Dry in a vacuum oven at  $65^{\circ}\text{C}$  for ~5 hours. Keep covered and store in a desiccator to cool.
- 6.4. Re-weigh the cooled sample and calculate the weight loss as the moisture.

## 7. Free Water Determination:

- 7.1. Weigh the empty aluminum dish and record the weight.
- 7.2. Mix sample thoroughly and weigh  $20.0\text{g} \pm 0.1\text{g}$  into the aluminum dish.
- 7.3. Dry at  $45^{\circ}\text{C}$  for a minimum of 8 hours or overnight. Keep covered and store in a desiccator to cool.
- 7.4. Re-weigh the cooled sample and calculate the weight loss as the free water.

## 8. Calculations:

For samples not dried prior to moisture / free water analysis:

$$\% \text{ Moisture} = \frac{(W - D) * 100}{(W - A)}$$

For samples pre-dried prior to moisture analysis

$$\% \text{ Moisture} = \frac{(100 - P) * (W - D)}{(W - A)} + P$$

Where:

- W = The weight of the wet sample and the aluminum dish
- D = The weight of the dried sample and the aluminum dish
- A = The weight of the aluminum dish (including sand, if used)
- P = Pre-dry % (for samples pre-dried prior to moisture analysis)

**Approvals:**

**Written By:**

Stacy Aylesworth  
Stacy Aylesworth  
Senior Environmental Scientist

2/15/21  
Date

**Approved By:**

Maryam Khosravifard  
Maryam Khosravifard  
Environmental Program Manager I

2/19/21  
Date

Sarva Balachandra  
Sarva Balachandra  
Quality Assurance Officer

3/10/21  
Date

