

Urea Nitrogen in Fertilizer Samples

1. Scope:

This procedure is to be used for the analysis of urea nitrogen in fertilizer.

2. Principle:

The amount of urea is determined by analyzing samples for ammoniacal + urea nitrogen. Ammoniacal nitrogen is determined separately and is subtracted from the ammoniacal & urea result. Samples are prepared as described in RA-SP-SMPL-PREP.

3. Safety:

The SDS for all chemicals shall be read before performing any part of this method. Gloves, safety glasses, and a lab coat shall be worn when handling hazardous materials and reagents.

Sulfuric acid is very toxic and extremely corrosive and shall be used in a fume hood. The acidification of samples containing reactive materials may result in the release of toxic gases and can be exothermic. Sample acidification and digestion shall be performed carefully in a fume hood.

4. Equipment, Reagents, and Supplies (equivalents are acceptable):

- 4.1 Kjeldahl distillation Unit
- 4.2 Kjeldahl flasks – 800 mL
- 4.3 Rubber stoppers – Size 7
- 4.4 Erlenmeyer flasks – 500 mL
- 4.5 Urease powder (Fisher Scientific cat# U21)
- 4.6 Magnesium oxide (Fisher Scientific cat# M68-3)
- 4.7 Calcium chloride solution prepared by adding 25g CaCl_2 in a 100mL volumetric flask and filling to the line with water
- 4.8 2% Boric acid solution with indicator (Fisher Scientific cat# 1064-1)
- 4.9 Defoaming solution (Spectrum cat# A1302)
- 4.10 Standardized sulfuric acid solution - 0.5N (Fisher Scientific cat # SA215)
- 4.11 Phosphate buffer prepared by adding 13.61g KH_2PO_4 and 17.42g K_2HPO_4 in a 200mL volumetric flask and filling to the line with water

5. Sample Analysis

- 5.1 For samples containing <5% urea, weigh 2g into a labeled Kjeldahl flask (record the weight). For samples containing $\geq 5\%$ urea, weigh 1g of sample.
- 5.2 Prepare a blank consisting of 250mL water and process using the same procedure as the samples.
- 5.3 Add 0.25g urease and 10mL phosphate buffer to the Kjeldahl flasks. Wash down the sides of each flask with 250-300mL water (except for blank), stopper tightly, and gently swirl to mix. Let stand ≥ 1 hour at room temperature.
- 5.4 For each sample, add 100mL 2% boric acid solution with indicator to a labeled 500mL Erlenmeyer flask and place beneath each unit of the Kjeldahl distillation unit with the glass tube immersed into the boric acid solution.
- 5.5 Turn on the cooling system for the Kjeldahl unit and turn on the burners for the units that will be used.
- 5.6 Add a few drops of defoaming solution (enough so foaming stops), 5mL CaCl_2 solution, and $\geq 2\text{g}$ MgO to the Kjeldahl flask. Rinse the neck with water and immediately connect the flask to the condenser of the Kjeldahl unit. Gently swirl the flask to mix the contents.
- 5.7 Distill $\sim 100\text{mL}$ (100mL minimum) into the boric acid solution. The total volume will be $\sim 200\text{mL}$.
- 5.8 Titrate with standardized sulfuric acid solution to the endpoint (color change to bright red).
- 5.9 The result includes both ammoniacal and urea nitrogen.

6. Calculations:

$$\% \text{ ammoniacal \& urea} = \frac{(S - B) * 0.5 * 1.4007}{\text{Sample Weight (g)}}$$

Where:

- S = Amount (mL) of sulfuric acid used to titrate sample
- B = Amount (mL) of sulfuric acid used to titrate blank
- 0.5 = Normality of the sulfuric acid solution
- 1.4007 = Conversion factor

Use the ammoniacal nitrogen method to determine the amount of ammoniacal nitrogen present in the sample.

$$\% \text{ urea} = \% \text{ ammoniacal \& urea} - \% \text{ ammoniacal}$$

7. References:

- 7.1 Official Methods of Analysis of AOAC, 17th ed., Method 941.04, Chapter 4.3.01
- 7.2 Official Methods of Analysis of AOAC, 19th ed., Method 955.04, Chapter 2.4.03
- 7.3 Standard Methods for the Examination of Water and Wastewater, 20th ed., Method 4500-NH₃ C, Chapter 4 – 105
- 7.4 Methods of Soil Analysis, Part 3 - Chemical Methods, SSSA, Regular Kjeldahl Method, page 1104

Approvals:

Written By:

Stacy Aylesworth
Stacy Aylesworth
Senior Environmental Scientist (Supervisor)

2/15/21
Date

Approved By:

Maryam Khosravifard
Maryam Khosravifard
Environmental Program Manager I

2/19/21
Date

Sarva Balachandra
Sarva Balachandra
Quality Assurance Officer

3/16/21
Date

