

CULTURAL PROCEDURES - GENERAL REQUIREMENTS

[Unless otherwise stated all tolerances are $\pm 5\%$]

1. Work Area

- a. Level table or bench, ample working space and utilities
- b. Clean, well ventilated, temperature 16-27C reasonably free from dust and drafts
- c. Well-lighted, > 50 foot-candles at working surface (pref. 100)
- d. Microbic density of air @ 15 colonies/plate in 15 min exposure, or @ 10 colonies/PAC plate in 15 min exposure, if not corrective actions taken
- e. Freedom from congestion and traffic, only compatible laboratory functions performed
- f. Safe working environment - Refer to OSHA
 - 1. Eating and drinking not permitted in laboratory
 - 2. Food and drinks for consumption not stored in laboratory
 - 3. Analysts wear buttoned/snapped lab coats/uniforms and protective eye-wear, lab coats/uniforms remain on-site
 - 4. Safety equipment available
 - 5. MSDS sheets in laboratory available to analysts
 - 6. Has functioning fume hood with acceptable sash (if necessary, see DMSCC procedure)
 - 7. Flammable solvent areas continuously well ventilated and temperature controlled
 - 8. Proper disposal of potentially hazardous materials
 - a. Contaminated samples disposed of properly
 - b. Contaminated glassware or plasticware disposed of or decontaminated properly
 - c. Hazardous chemical disposed of properly
- g. Storage Space
 - 1. Cabinets, drawers, and shelves adequate
- h. Areas neat, clean and orderly

- i. Floors clean, walls and ceilings in good repair _____
- j. Laboratory free of insects and rodents _____
- 2. Records _____
 - a. All laboratory related records maintained and available for announced surveys _____
 - 1. Three (3) years for state central labs _____
 - 2. Two (2) years for other labs, minimum requirement, States may require longer periods) _____
 - b. Quality control and sample records available to laboratory evaluation officer during survey _____
 - c. Records contain written corrective actions when taken _____
 - d. Records written in ink or other indelible substance, pencil or erasable ink not allowed _____
 - e. Corrections to quality control records, bench sheets and reports follow the requirements below: _____
 - 1. Make a single line through the incorrect information _____
 - 2. Write in the correct information next to the incorrect information _____
 - 3. Person making the correction initials the information _____
 - 4. If not obvious, include reason for correction _____
 - f. Requirements for electronic/computer records _____
 - 1. Software must be well documented _____
 - 2. Protocols and policies must be clearly documented _____
 - 3. Records must be indexed and cross referenced to allow easy review, or must be printed and made available _____
 - 4. Records must be secure from unauthorized access and changes _____
 - 5. When corrections are necessary the old information must be retained, the person making the correction must be identified and the reason for the change recorded _____
 - 6. If records are not available at time of audit, facility will be cited for not having records and will be subject to penalties _____

APPARATUS & MATERIALS

3. Thermometers

- a. National Institute of Standards and Testing (NIST)
Certified Thermometer, or equivalent, with
certificate Serial Number _____

- 1. Graduation interval not more than 0.5C (0-100C)
otherwise not more than 1.0C (< 0 or > 100C) _____
- 2. Calibration date on certificate _____
- 3. Annually, checked at the ice point Date _____
- b. Range of test thermometers appropriate for designated
use _____
- 1. Mercury-in-glass, alcohol/spirit or digital in
degrees centigrade _____
- 2. Plastic lamination recommended for mercury
thermometers _____
- c. Graduation interval not more than 0.5C (0-100C)
otherwise not more than 1.0C (< 0 or > 100C) _____
- d. Accuracy of test thermometers checked against certified
thermometer _____
- 1. Accurate to ± 1 C when checked at temperature(s)
of use _____
- 2. Results recorded and thermometers tagged _____
 - a. Tag includes identification/location, date
of check, calibration temperature and
correction factor(s) (read to within ± 0.5 C) _____
- e. All test thermometers accuracy checked before initial
use and annually, including autoclave maximum registering
and hot air oven thermometers _____
- f. Electronic thermometers checked before initial use and
annually as described above _____
- g. Automatic temperature recording instruments, if used,
compared weekly against an accurate thermometer,
results recorded _____
- h. Dial thermometers not used in the laboratory _____

4. Refrigeration (Sample _____) _____
(Reagent _____) _____
- a. Size adequate for workload _____
 - b. Maintains samples at 0-4.4C; if temperature out of range, record samples as not analyzed (NA) _____
 - c. Used for storage of milk or milk products, media and reagents only _____
 - 1. Not to be used to store food or drink _____
 - d. Record temperature (corrected) daily, in AM and PM, from two thermometers with bulbs immersed in liquid (in sealed containers) _____
 - e. Thermometers located on upper and lower shelves of use _____
5. Freezer (_____)
- a. Size adequate for workload _____
 - b. Maintains -15C or below _____
 - c. Used for storage of frozen milk products, controls, media and reagents only _____
 - 1. Not to be used to store food or drink _____
 - d. Record temperature (corrected) daily, in AM and PM, thermometer with bulb immersed in antifreeze liquid (in sealed containers) _____
6. Pipets (Glass _____ Plastic _____ Pipettor _____)
- a. Appropriate capacity _____
 - b. Must conform to APHA specifications _____
 - c. Graduations distinctly marked with contrasting color _____
 - d. Discard those with broken tips, scratches or other defects _____
 - e. Pipettors, calibrated, fixed volume or electronic only _____
 - 1. Calibrate with ten (10) consecutive weighings once every 6 months (using separate tip for each weighing), average of all 10 weighings must be $\pm 5\%$ of specified delivery volume (by weight, or ~~0.1~~ 1.0 mL by volume using class A graduated cylinder), records maintained _____

2. Or, calibrate with 10 consecutive readings once every 6 months using the Artel PCS Pipette Calibration System, average of all 10 readings must be $\pm 5\%$ of specified delivery volume, records/printouts maintained _____
- a. Instrument, printer connected by manufacturer's supplied cable or instrument connected to computer via serial cable _____
- b. Instrument and printer (if applicable) connected to 120v/60Hz power _____
- c. Reagent kits and Instrument Calibrator kits stored at room temperature _____
1. Lot # _____ Exp. Date _____ _____
- d. Reagent Blanks and Sample solutions are the same lot _____
- e. Certificates of Calibration for Reagent Kit and Instrument Calibrator kit maintained in records _____
- f. Instrument Validation Guide available _____
- g. PCS Pipette Calibration System Procedure, follow manufacturer's Procedure Guide and instrument prompts _____
1. Uncover and insert Blank into the instrument _____
2. Determine which volumes are to be calibrated _____
3. Select the correct Sample Solution and aliquot sufficient amount into working vessel provided _____
4. Using the Pipettor to be verified, aspirate the Sample Solution from the working vessel and deliver it into the Blank seated in the instrument _____
5. When appropriate number of data are collected, press 'End of Run' button _____
6. Record results and file Pipette Calibration Certificate (printout) _____
- h. PCS Pipette System Quality Control _____
1. Following manufacturer's Procedure Guide and instrument prompts, perform an instrument calibration every 30 days or just prior to use _____

2. Record results and file Calibration Certificate (printout) _____
- i. PCS Calibration System Validation _____
 1. Upon receipt, validate the instrument by following the manufacturer's protocol _____
 3. Pipettors etched with identification (imprinted serial numbers acceptable) and tagged with date calibrated _____
 4. Tips (sterile for plate counts) appropriate to pipettor(s) being used _____
7. Pipet Containers _____
 - a. Used for sterilization, storage; non-toxic _____
8. Dilution Bottles and Closures, reusable _____
 - a. Bottles of borosilicate glass _____ or approved plastic _____ with smooth tops _____
 - b. Capacity 150 mL, indelibly marked at 99 ± 1 mL level _____
 - c. Closure non-toxic rubber stopper or plastic screw cap with liner _____
 - d. New Bakelite type plastic caps and closures treated to remove toxic residues, tested using a B. stearothermophilus type assay _____
 - e. Discard bottles and caps with chips, cracks, scratches or other defects _____
9. Petri Dishes (Glass _____ or Plastic _____) _____
 - a. Bottom at least 80 mm I.D., and 12 mm deep for plate counts _____
Brand _____
 - b. Bottom 86.1 - 87.0 mm I.D., and 12 mm deep for BsDA _____
Brand _____
 - c. Bottom flat and free from bubbles, scratches, or other defects _____
10. Petri Dish Container _____
 - a. Used for sterilization, storage; non-toxic _____
11. Hot-Air Sterilizing Oven (_____) _____
 - a. Sufficient size to prevent crowding of interior in normal usage _____

- b. Constructed to provide uniform temperature in chamber _____
 - c. Thermometer or temperature recorder with adequate range (to 220C) _____
 - 1. Thermometer checked at temperature of use for accuracy before initial use, records maintained _____
 - 2. Thermometer bulb immersed in sand _____
 - d. Records maintained for each sterilization cycle including date, start-up time, time sterilization temperature reached, and length of time at sterilization temperature _____
 - e. Temperature indicator used each load _____
 - f. Performance checked with full load and recorded quarterly (preferably weekly) using spore (B. subtilus) strips, include positive control check, results maintained _____
 - 1. Brand: _____
 - 2. Lot #: _____ Exp. Date: _____
12. Sterilization by Dry Heat _____
- a. Material in center of load heated to ~~21~~ 170C for ~~2~~ hrs _____
 - b. Oven not crowded (< 75% of shelf in gravity type, 90% in forced air type) _____
13. Autoclave (Media _____) _____
- (Waste _____) _____
- a. Sufficient size to prevent crowding of chamber _____
 - b. Thermometer or temperature recorder-controller properly located to register chamber temperature _____
 - c. Has pressure gauge and properly adjusted safety valve _____
 - d. Connected to suitable saturated steam line or steam generator _____
 - e. Chamber temperature checked at least quarterly (preferably more frequently, ex. weekly with sterility check) with full load with maximum registering thermometer and results recorded _____
 - f. Cycle timing checked quarterly and found to be accurate, record maintained _____

- g. Records maintained for each sterilization cycle including date, start-up time, temperature and time temperature reached, length of time at temperature, time at end of run, time removed and item(s) autoclaved (including waste) _____
1. Strip recorders that provide the above information are acceptable if strips (or copies) are maintained in permanent record, include items autoclaved, time removed and initials _____
2. Circular charts must be interpreted and must have written records to verify the information stated above _____
- h. Temperature indicator used each load _____
- i. Performance checked with full load and results recorded weekly using spore (B. stearothermophilus) strips or suspensions, include positive control check, results maintained _____
1. Brand: _____
2. Lot #: _____ Exp. Date: _____
- j. Routine maintenance performed and records maintained _____
14. Sterilization by Moist Heat _____
- a. Media autoclaved at 120±1C _____
1. Dilution buffer blanks for 15 min (30 min optional) _____
2. Media for 15 min (sugar broths as per manufacturer instructions) _____
- b. Media autoclaved within 1 hr of preparation _____
- c. Dilution buffer autoclaved on same day prepared _____
- d. Stoppers or caps slightly loosened to permit passage of steam and air _____
- e. All air expelled from autoclave before pressure allowed to rise _____
- f. Autoclave will reach 120±1C within 15 min (5 min pref) of starting air-exhaust _____
- g. Properly operating and calibrated temperature gauge (not a pressure gauge) relied on to insure sterilization _____
- h. After sterilization, pressure gradually reduced (~ 15 min) and media removed promptly when atmospheric pressure is reached _____

- i. Total time in autoclave less than 1 hour _____
- 15. Incubator and/or Incubator Room (SPC, PAC and Coliform)
 - (#1 _____) _____
 - (#2 _____) _____
 - a. Sufficient size to prevent crowding of interior _____
 - b. Shelves placed to assure uniformity at $32C \pm 1C$ _____
 - c. Chamber temperatures measured by not less than two thermometers with bulbs immersed in liquid (in sealed containers) _____
 - d. Thermometer located on the top and bottom shelves of use _____
 - e. Temperature (corrected) recorded from each thermometer twice daily (AM and PM) _____
 - f. Agar (10 - 12 mL) in SPC plates and/or (1 mL) in PAC plates must not lose more than 15% weight after 48 hrs incubation _____
 - 1. Agar weight loss of SPC and/or PAC plates tested quarterly and results recorded _____
 - a. Test minimum of two (2) plates/films per shelf in use, one on each side of shelf, preferably test 10 plates evenly distributed throughout the incubator _____
 - 2. Corrective action taken when criteria not met and records of corrective actions maintained _____
 - a. If weight loss is out of compliance take corrective actions (humidify incubator, reduce air flow, etc.) and retest as above and record _____
 - b. Use more agar (15 - 20 mL), if this option used laboratory must document that this amount of agar is routinely used for plating _____
- 16. Colony Counter _____
 - a. Quebec dark-field model or equivalent with satisfactory grid plate _____
- 17. Hand Tally, accurate _____
- 18. pH Meter (Milk Lab _____) _____
 - (Media Prep _____) _____
 - a. Electronic only, readable to 0.1 pH units _____

- b. Daily calibration and slope records and maintenance
log maintained when in use _____
 - c. Record date electrodes (double junction reference pref)
put into service (write in QC record and tag probe) _____
19. pH Measurement _____
- a. All measurements made at room temperature _____
 - b. Instrument standardized with known buffer solutions _____
 - 1. Three commercially prepared standard solutions used _____
 - 2. Each aliquot used once and discarded _____
 - 3. pH 4, 7 and 10 suggested for linearity and proper
function of meter _____
 - 4. Slope determined (95 - 102%) _____ each
time meter calibrated, records maintained _____
 - c. Medium pH recorded each time measured _____
 - d. Final (after sterilization) pH of each batch of
medium determined before use, records maintained _____
 - 1. Standard Methods Agar, pH 7.0 ± 0.2 _____
 - 2. Violet Red Bile Agar, pH 7.4 ± 0.2 _____
 - 3. Brilliant Green Bile Broth, pH 7.2 ± 0.2 _____
 - 4. PM Indicator Agar, pH 7.8 ± 0.2 _____
 - 5. Buffered rinse solution, 7.2 ± 0.2 _____
 - 6. Nutrient broth, pH 6.8 ± 0.2 _____
 - 7. Letheen Broth, pH 7.0 ± 0.2 _____
 - 8. Lauryl Tryptose Broth (LST), pH 6.8 ± 0.2 _____
 - 9. M-Endo Agar or Broth, pH 7.2 ± 0.2 _____
 - 10. MMO-MUG Medium, pH 7.4 ± 0.2 _____
 - 11. Stock phosphate buffer, pH 7.2 ± 0.2 _____
 - 12. Dilution buffer, pH 7.2 ± 0.2 _____

20. Balance (Milk _____) _____
(Media _____) _____
(Analytical _____) _____
- a. Electronic only, sensitive to ~~0.1~~ 0.1g for general laboratory purposes and proper sensitivity for calibrations and antibiotics _____
- b. Class S or S1, or equivalent ASTM 1, 2, or 3, weights _____
1. Certificate or other verification of authenticity _____
2. Free from excessive wear, filth and corrosion _____
3. Weights within class tolerance _____
- c. Checked monthly with weights corresponding to normal use of balance (ex. 100, 200, 500, 1000 mg, etc. for analytical balances, and 5, 10, 25, 50, 100, 150g, etc. for other balances), records maintained _____
- d. Checked at least annually, or when weights out of tolerance, by a qualified representative for good working order with proof of check in laboratory _____
1. Date of last check _____
21. Water Baths _____
- a. Thermostatically controlled to appropriate temperature(s) _____
- b. Water circulation capability, baths up to 64C _____
- c. Appropriate size for work loads _____
- d. Suitable water level maintained _____
22. Mechanical Dilution Bottle Shaker _____
- a. Type described in SMEDP, 11th Edition _____
- b. Other acceptable _____
23. Microwave Oven for Melting Media _____
- a. Analysts instructed to take extreme caution as media expands rapidly at the boiling point _____
24. Microbiologically Suitable (MS) Water _____
- a. Type _____
- b. System used _____

c. Monthly testing criteria

1. Standard plate count < 1,000 colonies/mL(< 10,000 colonies/mL if stored) _____

2. Total chlorine residual negative, recorded as less than the detection limit of test used _____

3. Resistance exceeds 0.5 megohm/cm or conductance is less than 2.0 umhos/cm (at 25C) _____

a. Brand: _____ Std. _____

b. Test performed in another lab _____

d. Tested annually for total metals (Pb, Cd, Cr, Cu, Ni and Zn), not to exceed 0.05 mg/L for each metal and not to exceed 0.1 mg/L total for all metals _____

e. If criteria not met, corrective action(s) taken and recorded in QC record _____

f. Records maintained _____

25. Dilution Buffer and Blanks _____

a. Stock phosphate buffer (Date prepared _____) _____

1. Prepared in laboratory (34g KH_2PO_4 /liter) with MS water _____

2. Purchased commercially prepared _____

3. Lot No. _____ Exp. Date _____

Rcd. Date _____ Date Opened _____

4. Place in small containers (200 mL), autoclave and store in refrigerator _____

b. Stock MgCl_2 Solution, Optional (Date prepared _____) _____

1. Prepared in laboratory (38g MgCl_2 /liter or 81.1 g $\text{MgCl}_2 \cdot 6\text{H}_2\text{O}$ /liter) with MS water _____

2. Purchased commercially prepared _____

3. Lot No. _____ Exp Date _____

Rcd. Date _____ Date Opened _____

4. Place in small containers (200 mL), autoclave and store in refrigerator _____

c. Prepare dilution buffer with 1.25 mL stock buffer/liter of MS water _____

1. Optionally, add 5 mL of stock MgCl_2 /liter of MS water _____

- d. Dilution bottles filled to contain 99 ± 2 mL dilution buffer after sterilization _____
1. After sterilization and after cool visually observe and discard any blanks with < 97 or > 101 mL _____
 2. Of remaining blanks appearing to have the correct volume, check 1 blank for every 25 that were made using a class A graduate cylinder (or equivalent) _____
 3. Maintain records of volume checks, including batch size _____
 4. If any blanks out of tolerance, discard entire lot, record lot as discarded _____
- e. Blanks tested at 6 month intervals for toxic substances _____
1. Plate milk dilution at 0, 15, 30, 45 min _____
 2. If the 45 min count is 20% less than 0 min count, determine cause and retest after correction made, records maintained _____
- f. Alternatively, commercially prepared dilution buffer blanks used Brand _____
- Lot No. _____ Exp. date _____
- Rcd. Date _____
1. Volume records maintained as above _____
 2. Toxicity checked as above on each new lot received _____
 3. Check pH and record _____
- g. Records maintained _____
- h. Corrective action taken when criteria not met, records maintained _____
26. Reagent Chemicals - of ACS Grade _____
27. Media _____
- a. Use dehydrated medium of correct composition _____
1. Each bottle dated on receipt (in lab or by central receiving, which ever first) and when first opened for use _____
 2. Stored as specified by manufacturer; after opening, each bottle tightly capped following each use _____
 3. Commercially sealed medium kept no longer than manufacturer's expiration date _____

4. Opened bottles used until manufacturer's expiration date _____
5. Discarded if any change is noted in appearance or hydration regardless of manufacturer's expiration date _____

b. Plate Count Agar _____

1. Composition Pancreatic Digest of Casein 5 g
 Yeast Extract 2.5 g
 Glucose 1 g
 Agar 15 g
 MS water to make 1 L

2. Lot No. _____ Exp. Date _____
Rcd. Date _____ Date Opened _____

c. Petrifilm Aerobic Count (PAC) Plate _____

1. Lot No. _____ Exp. Date _____
Rcd. Date _____ Date Opened _____

d. Violet Red Bile Agar _____

1. Composition Yeast Extract 3 g
 Peptone or Gelysate 7 g
 Bile Salts 1.5 g
 Lactose 10 g
 Sodium Chloride 5 g
 Neutral Red 0.03 g
 Crystal Violet 0.002 g
 Agar 15 g
 MS water to make 1 L

2. Boil 2 min, temper and use within 3 hours (do not autoclave) _____

3. Lot No. _____ Exp. Date _____
Rcd. Date _____ Date Opened _____

e. Petrifilm Coliform Count (PCC) Plate _____

1. Lot No. _____ Exp. Date _____
Rcd. Date _____ Date Opened _____

f. Petrifilm High Sensitivity Coliform Count (HSCC) Plate _____

1. Lot No. _____ Exp. Date _____
Rcd. Date _____ Date Opened _____

g. Brilliant Green Lactose Bile Broth _____

1. Composition Peptone or Gelysate 10 g
 Lactose 10 g
 Oxgall 20 g
 Brilliant Green 0.0133 g
 MS water to make 1 L

2. Lot No. _____ Exp. Date _____
Rcd. Date _____ Date Opened _____

h. PM Indicator Agar _____

1. Composition

Beef Extract	3 g
Peptone	5 g
Tryptone	1.7 g
Soytone	0.3 g
Dextrose	5.25 g
Sodium Chloride	0.5 g
Dipotassium Phosphate	0.25 g
Polysorbate 80	1 g
Brom Cresol Purple	0.06 g
Agar	15 g
MS water to make	1 L

2. Lot No. _____ Exp. Date _____
Rcd. Date _____ Date Opened _____

i. Buffered Rinse Solution _____

1. Composition

Stock Phosphate Buffer	1.25 mL
10% Na Thiosulfate Solution	5 mL
Azolectin	4 g
Tween 20	10 g
MS water to make	1 L

2. Weigh hygroscopic Azolectin rapidly and dissolve by heating over boiling water _____

3. Date prepared _____

j. Nutrient Broth (laboratory use only) _____

1. Composition

Beef Extract	3 g
Peptone	5 g
MS water to make	1 L

2. Lot No. _____ Exp. Date _____
Rcd. Date _____ Date Opened _____

k. Letheen Broth _____
(For use with Petrifilm, Do not use diluents containing thiosulfate or sodium citrate)

1. Composition

Peptamin	10 g
Beef Extract	5 g
Lecithin	0.5 g
Sorbitan Monooleate	5 g
Sodium Chloride	5 g
MS water to make	1 L

2. Lot No. _____ Exp. Date _____
Rcd. Date _____ Date Opened _____

- o. MMO-MUG Medium _____
1. Commercial or lab prepared media containing MMO-MUG _____
2. Composition
- | | |
|---|----------|
| Ammonium Sulfate | 5 g |
| Manganese Sulfate | 0.0005 g |
| Zinc Sulfate | 0.0005 g |
| Magnesium Sulfate | 0.1 g |
| Sodium Chloride | 10 g |
| Calcium Chloride | 0.05 g |
| Sodium Sulfite | 0.04 g |
| Amphotercin B | 0.001 g |
| o-nitrophenyl- β -D-galactopyranoside | 0.5 g |
| 4-methylumbellifery- β -D-glucuronide | 0.075 g |
| Solanium | 0.5 g |
| Hepes Buffer | |
| Sodium Salt | 5.3 g |
| Organic Acid | 6.9 g |
| MS water to make | 1 L |
2. Lot No. _____ Exp. Date _____
 Rcd. Date _____ Date Opened _____
- p. Charm E*Colite _____
1. Lot No. _____ Exp. Date _____
 Rcd. Date _____
28. Medium Preparation _____
- a. Media-making utensils borosilicate glass, stainless steel, or other non-corrosive equipment _____
- b. Weigh required amount of dehydrated medium or ingredients _____
- c. Combined with required amount MS water, dissolved and mixed in a suitable container _____
- d. pH adjusted if necessary _____
- e. Heated (covered), not under pressure, if necessary, to complete solution (microwave preparation not allowed) _____
- f. Water restored, as necessary, to compensate for loss due to evaporation _____
- g. Distributed into suitable containers so that no part of medium is more than 2.5 cm from any surface _____
1. In general, containers filled no more than half of total volume _____
- h. Suitable container closures used and autoclaved as necessary _____

29. Prepared Media Storage _____
- a. Protected from water loss and light _____
 - b. Only screw-capped containers kept no more than 6 months _____
 - c. Prepared Charm PMI plates, kept no more than 5 days in sealed container at 0-4.4C (tag with date of preparation) _____
 - d. BGB broth at room temperature _____
 - 1. Screw capped tubes for 3 months _____
 - 2. Loose (slip) capped tubes for 1 week _____
 - 3. Stored in dark _____
 - e. Petrifilm plate storage _____
 - 1. Refrigerate unopened packages of Petrifilm plates at or below 8C, if frozen allow 30 min room temperature thaw time before opening packages _____
 - 2. Use before expiration date on package _____
 - 3. After opening, return unused plates to foil pouch, seal pouch by folding and taping/clipping open end shut _____
 - 4. Store re-sealed packages ~~at~~ 21C, ~~at~~ 50% relative humidity. **Do not refrigerate opened packages.** _____
 - 5. Use Petrifilm plates within one month after opening package (tag with date opened) _____
 - f. Pre-dispensed rinse solutions for containers _____
 - 1. Dispense in appropriate volume (20, 50, 100 mL, or other) and sterilize _____
 - 2. Perform quality control checks for volume (100±2 mL) as described in cultural procedures item 25d _____
30. Detergent Suitability Test _____
- a. Detergent residue test performed if laboratory washes and re-uses glassware (not required if only disposable items used) _____
 - b. Detergent is suitable for laboratory use _____
Brand _____ Brand _____
 - c. Test each new brand/lot, records maintained _____
31. Cleaning Pipets _____
- a. Used pipets discarded in disinfectant _____

- b. Rinsed in tap water at 15-30C _____
- c. Thoroughly washed with suitable detergent and rinsed _____
- d. Cleaned with strong cleaning solution such as acid dairy cleaner as necessary _____
- e. Final rinse with MS water _____
- f. Several pieces from each batch tested (preferably while still wet) for residual acid or alkali with aqueous 0.4% bromthymol blue. If color reaction not dark green to light blue, re-rinse and test again. Records maintained _____

32. Cleaning Other Glassware and Apparatus _____

- a. Heated to 85C or disinfected unless pathogens suspected; then sterilization required prior to washing _____
- b. Washed with hot water and suitable detergent and rinsed _____
- c. Machine washed (_____) _____
- d. Hand washed _____
- e. Final rinse with MS water _____
- f. Several pieces from each batch tested (preferably while still wet) for residual acid or alkali with aqueous 0.04% bromthymol blue. If color reaction not dark green to light blue, re-rinse and test again. Records maintained _____

SAMPLES

33. Laboratory Requirements _____

- a. Record time, date, and temperature of samples as received _____
- b. Determine sample temperature _____
 - 1. Insert a pre-cooled thermometer into temperature control (pre-cooling of electronic/digital thermometer probes is not necessary) _____
 - 2. Temperature control must be at least half the size of the largest test container _____
 - 3. Performed by trained personnel, not by collector, establish record to indicate training performed _____
- c. Do not accept samples if temperature control for each tank truck and each plant or delivery truck group of samples is missing _____
- d. Do not accept or test samples if sample containers are leaking _____

- e. Do not accept or test samples if samples are unprotected and/or submerged in ice/ice water slush _____
- f. Do not accept fluid samples, which are frozen, for microbial or somatic cell analysis _____
- g. Do not accept raw samples if sample containers have no head space (about $\frac{3}{4}$ full) _____
- h. If milk sample temperature control exceeds 4.4C on receipt, do not test microbiologically (samples may be tested if temperature does not exceed 7C **and** time of receipt is ~~03~~ hours from collection and sample receipt temperature is no greater than that at collection) _____
- i. Store fluid samples at 0-4.4C until tested, if storage temperature exceeds 4.4C prior to testing record as LA _____
- j. Do not proceed with analysis if above criteria have not been met _____
- k. Record date, time and temperature of samples when tested _____
- l. Begin testing of samples within 48 hr. of first collection (if time of collection not available use 12:01 am of date of collection) _____
- m. If chemical tests are made, remove portions for microbial analyses first _____

34. Sample Bench Sheet Requirements _____

- a. Must show date, time and temperature collected, along with name of official sampler _____
- b. Must show date, time and temperature when brought to the laboratory, along with whom received the samples _____
- c. Must show date and time of analysis, temperature of samples at start of analysis, and names of analyst(s) performing test(s) _____
- d. Sample bench sheets or records must contain all results (raw and calculated in their proper format) of tests performed and the results of all controls that apply to each test _____
 - 1. Plate count procedure controls include:
 - a. Microbic air density _____
 - b. Dilution buffer _____
 - c. Pipets _____
 - d. Agar _____
 - e. Temperature of agar at plating ($45 \pm 1C$) _____
 - 2. Results of inhibitor test(s) accompany all plate count and coliform results _____

3. Control results recorded for each inhibitor test performed _____

4. All above recorded on sample bench sheets _____

MISCELLANEOUS

35. Laboratory Practices _____

a. Personnel adequately trained and/or supervised _____

b. Satisfactory participation in annual split samples _____

c. Copies of current, applicable FDA 2400 series survey forms in laboratory _____

d. Copy of current edition Standard Methods in laboratory _____

e. Copy of AOAC Manual of Methods in laboratory if necessary _____

f. Copy of written Quality Assurance Plan, required for state central laboratories _____

g. Laboratory management has signed and returned the agreement to abide by the provisions of the NCIMS and the procedures for the Evaluation of Milk Laboratories _____

h. Laboratory evaluation officer conducted survey unobstructed by laboratory or facility personnel _____