

**PHOSPHATASE TEST - CHARM PASLITE - ALKALINE PHOSPHATASE TEST
USING CHARM II 6000/6600 AND LUMINOMETER/LUMINATOR**

[Unless otherwise stated all tolerances are $\pm 5\%$]

SAMPLES

1. Laboratory Requirements (see CP, items 33 & 34) _____

APPARATUS

2. CP, items 1 - 32 (as necessary) _____

- a. Unless otherwise stated, "shake vigorously" refers to standard microbiological mixing, i.e., 25 times in a one foot arc in seven seconds _____

3. Incubator Block for 13 x 100 mm test tubes or 2 mL microtubes _____

- a. Thermostatically-controlled at $35 \pm 1^\circ\text{C}$ _____
- b. Temperature checked by electronic display or by thermometer in small well in block or by liquid immersion, records maintained _____

4. Pipettors and Pipets _____

- a. Fixed volume or electronic, 100 μL _____
- b. Calibration checked as specified in CP item 6e, records maintained _____
- c. Disposable, 10 mL (ASTM) pipet with 0.1 mL graduations _____

5. Reagent Dispenser _____

- a. Fixed volume or electronic, 1.0 mL _____
- b. Calibration checked (CP item 6e) with 10 mL ASTM graduated cylinder, records maintained _____

6. Test tubes or microtubes _____

- a. Test tubes for Charm II 6600/Charm II 6000 systems, disposable borosilicate glass 13 x 100 mm, dirt and scratch free _____
- b. Microtubes - for Luminometer/Luminator, 2 mL screw cap _____

7. 6000/6600 or Luminometer/Luminator Analyzer

a. Operating instructions available

b. Definitions:

1. Fluid white milks - including skim through whole fat milk
2. Unflavored liquid dairy products - including half and half, buttermilk, creams (light, medium and whipping), and etc
3. Flavored liquid dairy products - Liquid products that can be accurately pipetted, containing flavor additives and/or thickening agents including flavored milks, and etc
4. Solid/semisolid dairy products - thick dairy products not able to be pipetted, solid and/or powdered additives, including cheese, yogurt, cottage cheese, ice cream mix, ice cream milk, whey, sour cream, and heavy cream - 36%, and etc

8. Water Bath, circulating, $34 \pm 1^\circ\text{C}$ and $63 \pm 1^\circ\text{C}$ (or $66 \pm 1^\circ\text{C}$ if fat > 10%), or 13 x 100 test tube dry well heater blocks acceptable (Confirmation procedure)

9. Centrifuge - Charm II Heraeus (3,400 RPM), minifuge, or equivalent (1,200 - 2,000 g)

10. Handling and storage

a. Kit contains Reagent AP, Stopping Solution and Alkaline Phosphatase Positive Control Tablet

Kit: Lot # _____ Rcd. Date _____ Exp Date _____

1. For solid/semisolid dairy products, Diluent AP

Diluent AP: Lot # _____ Exp. Date _____

b. Reagents stored at $0-4.4^\circ\text{C}$ until expiration date

c. Stopping Solution must be at $18-24^\circ\text{C}$ at time of use, may be stored at room temperature, expiration date marked 2 month from room temperature storage

d. Bottles labeled with receive and open dates

CONTROLS

11. Negative Control

- a. Product type. Prepare at least 50 mL for calibrator and positive control preparation
 - 1. Fluid white milk - heat a sample of product (highest fat content) to $95 \pm 1^\circ\text{C}$ for 1 minute with stirring
 - 2. Flavored liquid dairy products - heat a chocolate (highest fat content) to $95 \pm 1^\circ\text{C}$ for 1 minute with stirring
 - 3. Unflavored liquid dairy products - heat pasteurized light cream to $95 \pm 1^\circ\text{C}$ for 1 minute with stirring
 - 4. Solid/semisolid dairy products - mix or knead 5 g of product (highest fat content) with 20 mL Diluent AP until homogeneous and heat to $95 \pm 1^\circ\text{C}$ for 1 min with stirring
 - 5. Note, if using 13 x 100 test tube dry well heater block at 95°C it takes 10 minutes to heat product to 95°C for one minute, use temperature control
- b. Cool rapidly in an ice bath and hold at $0-4.4^\circ\text{C}$,
- c. Kept at $0-4.4^\circ\text{C}$, the Negative Control may be used for up to 48 hours
- d. If desired, distribute 1 mL quantities in small tubes (Milk only), seal and freeze in a non-frost-free freezer, or place in a styrofoam container and place in the center of a frost-free freezer for no more than 2 months, vials labeled with preparation and expiration dates

12. Positive Control (for daily checks)

- a. Reconstitute positive control (450 mU/L) with negative control, item 11, as indicated on label, or alternatively use 350 mU/L calibrator (item 13a2a)
- b. Shake vigorously and let settle 10 minutes at $0-4.4^\circ\text{C}$ for re-suspension

1. For solid/semisolid dairy products only,
add 1 mL of reconstituted positive control
with 3 mL of Negative control to complete
preparation of positive control
- c. Shake vigorously again and use for test _____
- d. Positive controls and calibrators held at 0-4.4C may
be used for 48 hours, milk controls may be frozen at
-15C or lower for up to 3 weeks, thaw in refrigerator
prior to use _____
- e. With 6600 and C2Soft, enter either the triplicate RLU
average of Positive Control or triplicate RLU average
of 350 mU/L calibrator as the pos avg and CP in
C2Soft configuration file. Refer to C2Soft manual _____

CALIBRATION

13. With each new kit lot # check calibration of analyzer _____

- a. Prepare 350mU/L, 175mU/L, 44mU/L (milk only), 88mU/L
(flavored and unflavored only) calibrators using
Negative Control, item 11 _____
1. Rehydrate a calibrator tablet with 100 uL water,
mix to disperse tablet, wait 1 minute and mix
again _____
2. Add the specified volume of Negative Control to
each dissolved calibrator tablet to make calibrators:
 - a. Add 2.5 mL to make 350mU/L
 - b. Add 5 mL to make 175 mU/L
 - c. Add 10 mL to make 88 mU/L
(flavored and unflavored only)
 - d. Add 20 mL to make 44 mU/L
(fluid white milk only)
- b. Calibrate instrument by testing each calibration
control (350, 175, 44 (or 88) mU/L) in triplicate _____

6600 with C2Soft Software

- c. For fluid white milks, unflavored or flavored liquid
dairy product on the 6600 system with C2Soft software,
follow the Standard Curve Calibration procedure _____
1. Program has a separate assay line for each product
group, fluid white milk, flavored and unflavored
liquid dairy product _____

2. In calibrate mode, enter low concentration (44 or 88 mU/L) value, followed by 3 replicate counts _____
 3. Enter medium concentration (175 mU/L) value, followed by 3 replicate counts _____
 4. Enter high concentration (350 mU/L) value, followed by 3 replicate counts _____
 5. Calibration successful will be prompted at end of the procedure _____
- d. For solid/semisolid dairy products using the 6600 system with C2Soft, follow instructions for positive average or control point setup _____
1. Count 3 replicates of 350 mU/L control _____
 2. Control point is equal to average of triplicate counts _____

Luminometer/Luminator system

- g. For fluid white milk, unflavored or flavored liquid dairy products, determine average value for each calibrator _____
1. Set up a separate channel and calibration for each product group, fluid white milk, flavored and flavored liquid dairy products _____
 2. Check calibration _____
 - a. Average negative control tested in triplicate. Average must be less than 5 (less than 15 for flavored dairy products) _____
 - b. Average 44 mU/L (or 88 mU/L unflavored and flavored liquid dairy products) calibrator, must be between 35 - 53 mU/L (45 - 110 mU/L unflavored and flavored liquid dairy products) _____
 - c. Average 175 mU/L positive control, must be 145 - 205 mU/L _____
 - d. Average 350 mU/L calibrator, must be 320 - 400 mU/L _____
 3. If conditions are not met, recalibrate according to Luminometer/Luminator calibration instructions _____

- h. For solid/semisolid dairy products verify control point of 350 mU/kg _____
- 1. Count 3 replicates of 350 mU/kg control _____
- 2. Average 350 mU/kg positive control, must test 350 ± 105 mU/kg _____
- 3. If conditions are not met, recalibrate according to Luminometer/Luminator calibration instructions _____

DAILY PERFORMANCE CHECKS

14. Daily and to verify calibration, test a Negative Control (item 11) and Positive Control (item 12), for at least one product _____

- a. Test negative control beginning from item 15a or with solid and semi-solid dairy products from 15a2 _____
- b. Verify Negative Control calibration _____
 - 1. Fluid white milk test, unflavored and flavored assay value < 5 mU/L with luminometer or < 44 mU/L (<88 mU/L flavored and unflavored) with 6600 and C2Soft _____
 - 2. Solid and semi-solid dairy products < 350 mU/kg or less than or equal to the control point _____
- c. Verify Positive Control calibration _____
 - 1. Positive Control (450 mU/L) rehydrated with fluid white milk, flavored and unflavored fluid dairy products, must be 300-585 mU/L or 350mU/L calibrator must be 247-453 mU/L _____
 - 2. Solid and semi-solid dairy products, within $\pm 30\%$ of 350 mU/kg or the control point _____

TEST PROCEDURE

15. Procedure _____

- a. Prepare sample _____
 - 1. For fluid white milks, unflavored and flavored, invert filled retail container 25 times, each inversion a full cycle down and up, shake or vortex negative control _____

2. For solid/semisolid dairy products (**not including controls, item 11b & 12**) add 1 part to 4 parts Diluent AP _____
 - a. Mix or knead until homogeneous _____
 - b. Centrifuge for 3 minutes _____
 - c. Use liquid extract in item 15c _____
- b. Dispense 100 μ L of Reagent AP into test tubes or microtubes _____
- c. Dispense 100 μ L of the prepared sample (item 15a) or mixed controls (items 11d & 12) just above the Reagent AP and immediately mix _____
 1. Use a new pipet tip for each sample, place pipet tip in sample or prepared control (no more than 1 cm), draw up and remove tip from sample/control, expel once to pre-wet tip _____
 2. Draw sample or control into pipet tip, touch off to side of container _____
 3. Holding pipet 90 ∇ to lab bench at eye level, dry exterior of tip (if necessary) by wiping from the pipet toward the tip, be careful not to touch end of tip _____
 4. Dispense 100 μ L sample directly above surface of Reagent AP (do not dispense down side of test tube or microtube) _____
 5. Depress plunger several times to completely expel sample _____
 6. Mix test tubes or microtubes with a back-and-forth motion for 10 seconds, use a vortex mixer if available _____
- d. Place the test tube/microtube in the 35 \pm 1C incubator for 3 minutes _____
- e. Within 10 seconds after incubation add 1 mL of room temperature (18-24C) Stopping Solution _____
- f. Remove test tubes/microtubes from incubator, cap and shake each vigorously or vortex for 10 seconds _____

- g. Place test tube/microtube into analyzer within 3 minutes, tubes held at room temperature (Note: stability of count may be stabilized by placing tubes/microtubes in a room temperature bath)

1. **6600 with C2Soft software**

- a. Select appropriate assay type
- b. Enter ID of sample and press enter
- c. Load sample in analyzer and press enter
- d. In 5 seconds RLU reading will be displayed, mU/L value will appear in results or pop-up window
- e. For solid/semisolid dairy products, sample RLU will be compared to control point

2. **Luminometer/Luminator**

- a. Select AP Assay (or customized channel)
- b. Press Start or Enter
- c. In 5 seconds mU/L reading will be displayed

- h. Counting of all test tubes/microtubes must be completed in 3 minutes
- i. Samples with ≥ 350 mU/L (category 4) (or for solid/semisolid dairy products, values greater than or equal to control point) of ALP activity are suspect positive containing about 0.1% (v/v) or more raw milk and must be tested for microbial, and reactivated phosphatase (items 16 & 17)

CONFIRMATION

16. Microbial Phosphatase

- a. Heat 1.0 mL of suspect sample at $63 \pm 1^\circ\text{C}$ for 30 minutes, stirring or mixing every 10 minutes
 1. For semisolid/solid dairy products dilute 1.0 g suspect sample with 4.0 mL diluent AP, mix or knead until homogeneous
 2. If fat content is $> 10\%$, heat at $66 \pm 1^\circ\text{C}$ for 30 minutes

- b. Cool sample rapidly to 0-4.4C in an ice bath _____
- c. Test positive and negative controls following item 15 _____
- d. Test heated sample and unheated sample (original sample) following item 15 (semisolid/solid dairy products begin at item 15b) _____
- e. Interpretation _____
 - 1. Controls test as specified in item 14 _____
 - 2. If heated and unheated sample have equal activity ($\pm 30\%$, mU/L or RLU) the sample is regarded **Not Found** for residual phosphatase, the activity originally measured is microbial _____
 - 3. If the heated sample is more than 30% below unheated sample (mU/L or RLU), the sample contains milk phosphatase activity, either residual or reactivated _____

17. Reactivated Phosphatase _____

- a. Magnesium acetate solution commercially available _____
- b. Or, prepared in laboratory _____
 - 1. Dissolve 35.4 g of magnesium acetate tetrahydrate, $\text{Mg}(\text{C}_2\text{H}_3\text{O}_2)_2 \cdot 4\text{H}_2\text{O}$ in 25 mL MS water, warming slightly to aid dissolution _____
 - 2. Pour solution into 100 mL volumetric flask, rinse original container several times and add rinses to flask _____
 - 3. After cooling to room temperature, make up to 100 mL (stable for 1 year at 0-4.4C) _____
- c. Procedure _____
 - 1. Label separate test tubes as "Blank" and "Test" _____
 - 2. Add a 5.0 mL aliquot of sample (unheated, original sample not prepared as in 15a) to each test tube _____
 - a. For semisolid/solid dairy products, combine 2.5 g product and 10.0 mL Diluent AP _____
 - b. Mix or knead until homogeneous, and add 5.0 mL to clean test tubes labeled "Blank" and "Test" _____

3. Add 0.1 mL MS water to the sample labeled "Blank", and 0.1 mL magnesium acetate solution to the sample labeled "Test" _____
4. Cap tubes and heat both aliquots for 1 hr at 34±1C _____
5. Remove samples from water bath and cool rapidly to 0-4.4C in an ice bath _____
6. Dilute 1 mL of sample containing magnesium acetate (Test) with 5 mL (1:6 dilution) of negative control product (item 11), label tube as "Diluted Test" _____
7. Test undiluted sample containing no magnesium acetate (Blank) and diluted sample containing magnesium acetate (Diluted Test) for phosphatase activity following item 15 (semisolid/solid dairy products begin at item 15b) _____

d. Interpretation _____

1. If the diluted aliquot containing magnesium acetate (Diluted Test) has equal (±30%) or greater phosphatase activity than the undiluted aliquot containing no magnesium (Blank), the sample is regarded as **Not Found** for residual phosphatase, and the phosphatase originally measured is of reactivated origin _____

Dil. w/Mg (Test) - Undil. (Blank) = Reactivated

2. If the diluted aliquot (Diluted Test) contains less (30% below or less) activity than the undiluted aliquot (Blank) the sample is considered **Positive** for residual phosphatase _____

Dil. w/Mg (Test) < Undil. (Blank) = Residual

3. A false-positive for residual phosphatase may also be obtained if a reactivatable sample has been allowed to stand at elevated temperatures (20C) for periods of 1 hr or more before testing (SPC < 20,000/mL) _____

REPORT

19. Report as:

1. Residual phosphatase **Not Found (NF)**

- a. Report as < 44 mU/L (< 88 mU/L unflavored and flavored liquid dairy products or < 350 mU/kg solid/semisolid dairy products)

2. Residual phosphatase **Positive**

- a. Microbial and reactivatable phosphatase are not demonstrated
- b. Suspect positives greater than or equal ($\pm 30\%$) to 350 mU/L, category 4 or greater than the control point must be tested for microbial and reactivated phosphatase (items 16 and 17)
- c. Report mU/L values, mU/L range or greater than control point when equal to or greater than 44 mU/L fluid white milk, 88 mU/L unflavored and flavored liquid dairy products or 350 mU/kg semisolid/solid dairy products

3. Report as **Not Found** for residual phosphatase if:

- a. If microbial phosphatase present
- b. If reactivated phosphatase present
- c. If there is documentation to show that the product was treated such that reactivated phosphatase may be present